Spares

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bedfont
scientific contributions to health

Page

Introduction	1
Quick Start Guide	2
Pack Contents List	3
Instrument Layout	4
Maintenance/Troubleshooting	5
Operation	6
Calibration	10
Syringe Calibration	11
Specification	13
Returns Procedure	14
Spares & Warranty	15

Hydrogen (H_2) is generated in the intestinal lumen by bacterial action on carbohydrates in the large or small intestine. This resultant H_2 diffuses into the bloodstream and then to the alveoli, after which it can be detected in expiratory air. Levitt (1969)¹ demonstrated the correlation between intestinal lumen H_2 production and H_2 excretion in expiratory air. Thus, accurate measurement of H_2 in parts per million (ppm) in expiratory air reveals abnormal breakdown and/or malabsorption of carbohydrates.

The Bedfont Gastrolyzer provides a simple, user-friendly sampling technique with specific H_2 detection, which is reliable, inexpensive, and rapid².

Operation is straightforward. A T-piece sampling system enables end-expired breath to be sampled easily and hygienically, using one-time disposable cardboard tube mouthpieces. A modified low volume sampling technique enables children or people unable to give a sample to be monitored with a minimum of fuss, while neonates in incubators can be monitored simply by linking the Gastrolyzer to the expiratory limb of a ventilator. Alternatively, a face mask arrangement can be adopted.

The built-in 2-line LCD display shows the results within seconds, and stores the maximum value until 'switch off'. Thus technical analysis employing skilled personnel becomes redundant and, being fully portable with its own power source, the unit can be used in any location.

The Bedfont Gastrolyzer is supplied with full operating instructions, while maintenance amounts to no more than changing batteries when required, calibration at regular intervals and keeping clean by wiping down with a cleansing wipe.

The Gastrolyzer Portable Breath Hydrogen Monitor uses as its principle of detection a sealed electrochemical sensor which is specific to H_2 . The sensor is of the micro fuel cell type, designed to be maintenance free and stable over long periods of time. Because of the unique diffusion barrier, the sensor is unaffected by normal temperatures and has a linear response to H_2 concentration independent of pressure.

The Gastrolyzer can be used for the analysis of Carbohydrate Breakdown Deficiency Detection, Carbohydrate Malabsorption Detection, Lactose Intolerance Diagnosis, Bacterial Overgrowth Diagnosis and Determination of Time Passage Through Gut.

¹ Levitt, M.D. (1969): Production and excretion of Hydrogen gas in man. New Engl. J.Med. 281:122-127

² Fleming, S.C. (1990): Evaluation of hand-held Hydrogen monitor in the diagnosis of intestinal lactose deficiency. Ann.Clin.Biochem. 1990;499-500

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Voltooid aan de basiseisen van de Medical Device Directive
93/42/EEC Annex V. Certificate No. CE:01469.

1. Switch the unit ON, display will show:

**Bedfont
Gastrolyzer 2**

2. Unit will automatically enter a Zeroing mode, taking up to 2 minutes.
3. If Zeroing has been successful, the unit will beep and display will show:


**Ready
Gastrolyzer 2**

4. Connect T-piece and mouthpiece, ask patient to hold their breath for 15 seconds, and press the blue GO button.
5. Display now shows a 15-second countdown, with the last 3 seconds accompanied by beeps.
6. At the end of the countdown, the patient should exhale slowly, gently and deeply into the mouthpiece, and the display will show the ppm hydrogen reading.
7. The test is completed when the top line of the display shows "Ready" with the ppm reading on the bottom line:

**Ready
XXppm H₂**

8. Either switch the unit OFF, or press the blue GO button to start another test.
9. A new mouthpiece should be used for each person being tested.
10. The T-piece assembly should be removed between tests and the air around the sample vent stirred to ensure none of the previous sample remains in the sensor.

1. Gastrolyzer instrument
2. T-piece sampling system
3. Cardboard disposable mouthpieces (x3)
4. Calibration adaptor
5. Adult face mask
6. Child face mask
7. Infant face mask
8. Face mask sampling system
9. Calibration screwdriver
10. Battery
11. Operating manual
12. Carrying case

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
The 2-line message display provides user messages, warnings and displays the patient's end-tidal breath H₂ reading in ppm. In addition, the display shows a 15 second countdown during which the patient holds a lungful of breath before providing the end-tidal breath sample.

The Auto Zero enables the operator to zero the instrument before use, ensuring accurate readings every time.

Key

1. 2-Line LCD message display
2. T-piece sampling system
3. Disposable cardboard mouthpiece
4. ON/OFF Switch
5. ZERO button
6. GO button
7. Battery Compartment
8. Span adjustment
9. Calibration adaptor
10. Face mask sampling system
11. One-way valve mouthpiece

Routine Maintenance

- Calibrate the Gastrolyzer using Bedfont 200 ppm Hydrogen in air calibration gas every month (see Page 10). Once calibration has been carried out, the date should be recorded and the next calibration date diarised.
- Replace batteries when indicated.
- Replace T-piece sampling system every month or if visibly soiled or contaminated. It cannot be cleaned or sterilised. Record every time this is done and diarise the replacement date.

**Cleaning**

Wipe the Gastrolyzer instrument and external T-piece surfaces with a product specifically developed for this purpose. Bedfont provides an '**Instrument Cleansing Wipe**' (Part No.: **WIPE-V**)

NEVER use alcohol or cleaning agents containing alcohol or other organic solvents as these vapours will damage the Hydrogen sensor inside. Under no circumstances should the instrument be immersed in liquid or splashed with liquid.

Batteries

Batteries should be removed if the instrument is not likely to be used for some time.

Additional technical information can be made available on request; please contact Bedfont or its distributor.

Troubleshooting

There is no LCD or panel response under gas when sampling with new batteries?
Return the equipment immediately to Bedfont Scientific Ltd or the local representative. Under no circumstances attempt to locate and correct the fault yourself.

LCD readings are erratic and unsteady?

Check battery strength by switching the unit off and then back on again. If the LOW BATTERY warning message appears, switch unit off and replace battery.

The Gastrolyzer returns consistently low readings or none at all?

Check the condition of the sampling system. Under frequent use, saliva deposits can cause valves to stick. DO NOT attempt to clean them but replace the complete sampling system immediately. If, after replacement, there is no improvement in performance, return the Gastrolyzer to Bedfont Scientific Ltd or the local representative.

The LOW BATTERY warning message appears on the display?

Switch unit off and change battery immediately. Replacement battery must an alkaline PP3 9V battery or equivalent.

Your Gastrolyzer is dropped or immersed in liquid?

Return the instrument to Bedfont Scientific Ltd or the local representative for checking.

All accessories supplied are disposable. Replacement parts can be purchased from Bedfont Scientific Ltd or the local representative.

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Power On

Ensure battery is located in battery compartment. Switch ON/OFF switch to ON position. If low battery warning message is displayed, turn the unit off and replace the 9-volt alkaline battery, which is located in the battery compartment.

NOTE: The batteries should be removed if the instrument is not used for periods greater than one month.

After power up the display will show:

Bedfont
Gastrolizer 2

The top line then displays:

* ZEROING *

If appropriate the display may show:

SENSOR HIGH or SENSOR LOW

When Zeroing is complete the display will show:

SENSOR O.K.

This is followed by a long beep and the following display:

Ready
Gastrolizer 2

This zeroing process can take up to 2 minutes.

Manual Zeroing

The instrument can be manually zeroed at any time. Ensure that no gas is connected and the T-piece is removed before continuing. To perform a manual zero press the red Zero button. The top line of the display shows:

PLEASE WAIT

The bottom line shows:

Remove T-Piece

Followed by:

Ensure No Gas

Followed by:

* ZEROING *

When complete the instrument produces a short beep and displays:

Ready
Gastrolyzer 2

Taking a Test Reading

Connect the T-Piece sampling system with cardboard disposable mouthpiece attached. Ensure all connections and sampling system are pushed firmly in place to prevent leakage of sample.

Press the blue GO button at the same time asking the patient to hold their breath during the 15 second countdown. The display shows:

Breath Hold
XX secs

Where XX counts down from 15 to 0 in 1-second intervals. The last 3 secs are accompanied by short beeps and a final long beep. The top line of the display then shows:

!! Exhale !!

The patient should then exhale slowly and gently into the mouthpiece, emptying the lungs as far as possible.

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A powerful exhalation may dislodge the valve and will not necessarily provide the desired end-tidal air sample.

Alternatively if children or infants are to be monitored, a modified sampling technique can be used. This consists of a specially constructed breathing system, which allows exhaled breath to be directed to the instrument's sensor for analysis. The appropriate size face mask should be placed over the nose and mouth of the patient, care being taken to avoid leaks around the mask cushion. The patient can then breathe normally, whilst being encouraged to exhale as much as possible.

To use the face mask select the appropriate mask and attach to the face mask sampling system. This is then connected to the T-piece in place of the mouthpiece. See photograph opposite.

This procedure can start either after the 15-second countdown or during it. The display can then be observed to ensure the maximum reading is obtained.

NOTE: People with lung diseases or chest ailments may not be able to achieve the 15-second breath-hold. In such cases, initiate the 15-second countdown, wait for 0ppm to show on the LCD, and then instruct the patient to inhale. Hold breath as long as possible, and exhale through the Sampling System mouthpiece. Alternatively, if people with contagious diseases are being analysed, Bedfont recommend the Sampling System be replaced after use.

The top line of the display will show a bar graph of the real time reading and the lower line shows:

XX ppm H₂

Where XX is the peak held reading in ppm. This reading will be held until the instrument is switched off or another test is started.

When the test is completed the top display shows:

**Ready
XX ppm H₂**

Where XX is the peak held reading in ppm from the latest test. The instrument is now ready to take another test reading.

IMPORTANT: Always ensure the T-piece is removed between readings.

When starting a new reading, if the instrument detects the presence of gas in the sampling system the test will NOT be allowed to start, the instrument emits a low beep and the display will show:

!! No Zero !!
Ensure No Gas

If this occurs remove the T-piece and ensure the sensor is clear of gas before attempting another reading.

Taking a Reading with a Syringe (Low Volume Sampling Kit)

Before proceeding with this test method the instrument must be calibrated using a syringe see Syringe Calibration section

It is also possible to supply a sample using a syringe kit instead of directly from exhaled breath. Under this situation connect the sampling plug into the sensor housing of the instrument. Ensure it is firmly pushed in place to prevent leakage of the gas.

Connect the syringe to the sampling plug using the plastic tubing (see photograph opposite). Start the test in the normal way and when the display shows:

!! Exhale !!

Depress the syringe (note: you will feel a high back pressure, which is to protect the sensor) and maintain steady pressure. The instrument will then display the reading in the normal manner.

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For continued accurate use, calibration on a monthly basis is recommended. To calibrate the instrument a Bedfont calibration kit is required consisting of a gas can, fine control valve and flow indicator and plastic tubing.

Ensuring that the valve is in the OFF position (turn the valve knob clockwise), screw the fine control valve and flow indicator assembly to the gas can. This is best done by screwing the gas can into the valve.


With the tubing, connect the calibration adaptor and the flow indicator. Warming the end of the tubing using a hair dryer or lighter will assist connection. Insert the calibration adaptor into the T-piece sampling system in place of the mouthpiece. Insert the T-piece sampling system into the Gastrolyzer sensor housing, ensuring that all connections and sampling system are firmly pushed in place to prevent leakage of the calibration gas (see photograph opposite).

Open the fine control valve and allow the gas to flow at 0.5 litres/minute. To maintain this, adjust the flow so the ball in the Flow Indicator remains at the lower line (see diagram opposite).

To start the calibration process press the blue GO button followed by the red Zero button within 3 seconds. The display shows:

Calibrate
XXppm H₂

Where XX is the real-time ppm reading of the applied calibration gas, not the peak held value as with a normal reading.

Allow gas to flow through the instrument for 1½ minutes to ensure accurate calibration, again monitoring the rate of flow.

If after 1½ minutes, the ppm reading does not show between 195 and 205ppm, using the screwdriver, adjust the SPAN control on the underside of the instrument, until a reading of between 195 and 205ppm is displayed. Turning clockwise will increase the reading and turning anti-clockwise will decrease it.

When complete press the GO button to exit the Calibration mode.

Turn off the gas flow, remove T-piece sampling system and disconnect the calibration adaptor. Unscrew the fine control valve and flow indicator from the gas can and store safely. If valve is left in the can, the gas could escape.

If the normal method of taking readings is to use a syringe then it is recommended to calibrate the instrument using a syringe to ensure maximum accuracy.


The syringe calibration method is similar to the normal calibration method except that a syringe kit is required in addition to the normal calibration kit.

The syringe kit consists of a syringe, 2 pieces of tubing, a tubing T-piece and a sampling plug.

Connect the syringe to the calibration gas cylinder as indicated in the photograph opposite.

Connect the sampling plug into the sensor housing of the instrument. Ensure it is firmly pushed in place to prevent leakage of the gas.

Open the fine control valve and allow the gas to flow at 1.0 litres/minute. To maintain this, adjust the flow so the ball in the flow indicator remains at the centre line (see diagram below).

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
Draw a sample slowly into the syringe. Disconnect syringe from the T-piece sampling system and connect to sampling plug (see photograph below).

Turn off the gas flow, remove T-piece sampling system and disconnect the calibration adaptor. Unscrew the fine control valve and flow indicator from the gas can and store safely. If valve is left in the can, the gas could escape.

To start the calibration process press the blue GO button followed by the red Zero button within 3 seconds. The display shows:

Calibrate
XXppm H₂

Where XX is the real-time ppm reading of the applied calibration gas, not the peak held value as with a normal reading.

Depress the syringe (note: you will feel a high back pressure, which is to protect the sensor) and maintain steady pressure. If after the syringe has been emptied the ppm reading does not show between 195 and 205ppm, using the screwdriver, adjust the SPAN control on the underside of the instrument, until a reading of between 195 and 205ppm is displayed. Turning clockwise will increase the reading and turning anti-clockwise will decrease it.

When complete press the GO button to exit the Calibration mode.

Gas Detected:	Hydrogen (H_2)
Concentration Range:	0-250 ppm
Calibration Flowrate:	0.5 litres/minute (1.0 litre/minute with syringe)
Detection Principle:	Electrochemical sensor
Accuracy:	+/-2% of reading
Carbon Monoxide Cross Interference:	<15%
Alcohol Cross Interference:	Negligible
Display:	Liquid Crystal Display
Power Input:	PP3 alkaline battery (9v, 6LF22/Equiv)
Note: to avoid long start up times, ensure Gastrolyzer is switched OFF if battery is removed.	
Warm up Time:	Typically less than 120 seconds
Response Time:	Under 30 seconds to 90% FSD
Drift:	Typically less than 2% signal loss/month
Operating Temperature Range:	0-40°C
Operating Humidity Range:	0-100%
Sensor Operating Life:	2-3 years. 6 months warranty
Sensor Sensitivity:	1ppm
Dimensions:	63 (D) x 85 (W) x 144 (H)mm
Weight Including Battery:	Approx. 225 g
Construction:	Plastic
Recommended Storage Temp:	0-30°C

Meets the essential requirements of the Medical Device Directive 93/42/EEC Annex V. Certificate No. CE:01469.

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¹ Levitt, M.D. (1969): Production and excretion of Hydrogen gas in man. *New Engl. J.Med.* 281:122-127

² Fleming, S.C. (1990): Evaluation of hand-held Hydrogen monitor in the diagnosis of intestinal lactose deficiency. *Ann.Clin.Biochem.* 1990;499-500

- Please contact Bedfont Scientific Ltd before you return any monitors to us.
- You will need to have all the information about the monitor, including the serial number, and the fault, to give to the office.
- Bedfont will then issue you with a Returns Number.
- Please state this number on a compliment slip or headed paper when returning the monitor, ensuring your full details including telephone and fax numbers are clearly stated.
- We advise that you use a courier service when returning monitors, as this enables you to insure goods for loss or damage.
- When we have received your monitor we will send you an Order Acknowledgement.
- When the monitor has been examined we will send you an Engineer's Report, a quote for the repair, and an Authorisation Form.
- If your monitor is still in warranty we will repair it and return it to you with an Engineer's Report, free of charge. If it is found to simply need calibrating, you may be charged £40/€64/\$74 for the service, and the following procedure will apply.
- If you wish to proceed with the repair you will need to complete the Authorisation Form by ticking the relevant box for the repair and supply us with an **Official Purchase Order Number**, and fax it back to +44 (0)8700 844051.
- If you choose not to go ahead with the repair, there will be a handling fee of £30/€48/\$55.50. You will need to tick the relevant box on the Authorisation Form and return it to us. You must still supply us with an Official Order Number.
- As soon as Bedfont has received all the correct paperwork the unit will be returned to you. This will be at your cost if the unit is no longer in warranty.

Spares

<i>Part Number</i>	<i>Description</i>
EC60-T/P-V	Pack of 10 disposable T-pieces
EC60-MP-200	Box of 200 disposable mouthpieces
ONE-WAY MP/100-V	Box of 100 one-way valve mouthpieces
EC60-Adapt-V	Adaptor to allow new T-piece to be used with Gastrolyzer 1
EC60SSS	Syringe sampling system for low volume testing. Includes low volume plug, syringe and connectors.
EC60-ISS	Infant/child sampling system
EC60-IM	Infant face mask
EC60-MM	Medium sized face mask
EC60-AM	Adult face mask
EC60-CA	Calibration adaptor
EC60-COS	Replacement sensor
PP3-Batts	Replacement PP3 9v battery
EC60-Man	Replacement manual
EC-60-CC5	Carry case to house instrument and accessories
WIPE-V	Instrument Cleansing Wipes

The above spares are available from Bedfont Scientific Ltd, UK. For spares availability in all other countries contact your local distributor.

Warranty

Bedfont Scientific Limited warrants the Gastrolyzer 2 (batteries excepted) to be free of defects in materials and workmanship for a period of one year from the date of shipment. Bedfont's sole obligation under this warranty is limited to repairing or replacing, at its choice, any item covered under this warranty when such an item is returned intact, prepaid, to Bedfont Scientific Limited or the local representative.

Note: Sensors are guaranteed for a period of six months from the date of shipment from Bedfont.

These warranties are automatically invalidated if the products are repaired, altered or otherwise tampered with by unauthorised personnel, or have been subject to misuse, neglect or accident.

At the end of the product's life, do not dispose of any electronic instrument in the domestic waste, but contact Bedfont or its distributor for disposal instructions.

Bedfont Scientific Ltd

105 Laker Road,
Rochester Airport Industrial Estate
Rochester, Kent ME1 3QX England
Tel: +44(0) 8700 844 050
Fax: +44(0) 8700 844 051
E-mail: ask@bedfont.com
www.bedfont.com
www.gastrolyzer.com

bedfont

scientific contributions to health

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- Bitten setzen Sie sich mit Bedfont Scientific Ltd in Verbindung, bevor Sie ein Testgerät an uns zurücksenden.