



# NObreath®

FeNO Testing without limits.



FeNO monitoring  
made easy!

## Benefits of monitoring FeNO with the NObreath® :

- Non-invasive, quick and easy to perform<sup>1</sup>
- Aids in asthma management, assisting the correct prescription and making monitored adjustments
- Shows patient adherence to treatment<sup>4</sup>
- Aids in identifying patients who do/do not require on-going treatment<sup>2</sup>
- Aids in differentiating between allergic (eosinophilic) and non-allergic (non-eosinophilic) asthma<sup>3</sup>
- Shown to be superior to the majority of conventional tests of lung function, such as peak flow recording and spirometry<sup>1</sup>




Exclusive  
NObreath® forum




FREE FeNOchart™  
patient management  
software



Adult, child & ambient  
test modes



Create & save  
patient details



Onscreen animated  
flow meter for  
motivation

## Ideal for:

- GP's
- Respiratory Nurses
- Clinicians
- Medical Students

## Features and Benefits

Clinically proven



Low maintenance



Unlimited Usage\*



Adult & child testing



Easy to use



Low cost



## References

1. Andrew D. Smith, Jan O. Cowan, Sue Filsell, Chris MacLachlan, Gabrielle Monti-Sheehan, Pamela Jackson and D. Robin Taylor. Diagnosing Asthma: Comparisons between Exhaled Nitric Oxide Measurements and Conventional Tests. Am J Respir Crit Care Med Vol 169. pp 473-478, 2004.
2. D R Taylor, MW Pinenburg, A D Smith and J C D Jongste. Exhaled nitric oxide measurements: clinical application and interpretation. Thorax 2006;61:817-827.
3. Coumou HBel E. Improving the diagnosis of eosinophilic asthma [Internet]. Taylor and Francis online. 2017 [cited 15 March 2017]. Available from: <http://www.tandfonline.com/doi/full/10.1080/17476348.2017.1236688>
4. Beck-Ripp J, Griese M, Arenz S, Koring C, Pasqualoni B, Bufler P. Changes of exhaled nitric oxide during steroid treatment of childhood asthma. Eur Respir J 2002;19:1015-1019.

 [www.covita.net](http://www.covita.net) | 800.707.5751 

# Technical Specification

<b>Concentration range</b>		5-500ppb
<b>Display</b>		Full colour touchscreen
<b>Detection principle</b>		Electrochemical sensor
<b>Repeatability</b>		±5ppb of measured value ≤ 50ppb ±10% of measured value > 50ppb
<b>Accuracy</b>		±5ppb of measured value ≤ 50ppb ±10% of measured value > 50ppb
<b>Power</b>	<b>NObreath® monitor</b>	1 x main rechargeable Li-ion battery– Approx. 100 uses on fully charged battery Model:RRc1120. Voltage: 3.6V / 3.7V Capacity: 2350mAh/2000mAh 2x Li-ioncoin cell battery– Approx. 5 years Model: LIR2032 Voltage: 3.6V Capacity: 45mAh Model: LIR2450 Voltage: 3.7V Capacity: 120mAh
	<b>NObreath® Dock</b>	Mains powered Input: 5V, 0.5A Output: 5V, 0.5A
	<b>Plug</b>	Input: 100-240V ~ 50/60Hz., 0.2A Output: 5.0V, 1.0A
<b>T<sub>90</sub> response time</b>		≤10 seconds
<b>Temperature</b>	<b>Operating</b>	15-30°C (59-86°F)
	<b>Storage/transport</b>	0-50°C (32-122°F)
	<b>Calibration</b>	21°C ±4°C (17°C-24°C) 70°F ± 39°F (63°F - 77°F)
<b>Humidity</b>	<b>Operating</b>	20-80% RH (non-condensing)
	<b>Storage/transport</b>	5-95% RH (non-condensing)
<b>Operating/storage/transport Altitude</b>		-1700 ft.to 6300 ft.
<b>Sensor operating life</b>		5 years (subject to servicing)
<b>Limit of Detection</b>		5ppb
<b>Sensor drift</b>		<5% per annum
<b>Dimensions</b>		Approx. 90 x 159 x 59 mm (3.5 x 6.3 x 2.3 in)
<b>Weight</b>		Approx. 400g (0.9lb)
<b>Materials</b>	<b>NObreath® monitor</b>	Case: polycarbonate/ABS blend
	<b>NObreath® Dock</b>	SteriTouch® anti-microbial additive
<b>Breath test time</b>	<b>Adult</b>	12 seconds
	<b>Child</b>	10 seconds
	<b>Ambient</b>	30 seconds
<b>Warm-up time</b>		≤60 seconds
<b>Maximum ambient operating level</b>		350 ppb NO
<b>CO cross interference</b>		45ppm ≤17.6 ppb

\* Subject to correct use, maintenance and service. Tested up to 29,000 tests.

Visit [www.bedfont.com/resources](http://www.bedfont.com/resources) to view this document in other languages.



**Bedfont Scientific Ltd.**  
Station Road, Harrietsham, Maidstone  
Kent, ME17 1JA, England.  
Tel: +44 (0)1622 851122, Fax: +44 (0)1622 854860  
Email: [ask@bedfont.com](mailto:ask@bedfont.com) Web: [www.bedfont.com](http://www.bedfont.com)

© Bedfont® Scientific Limited 2022

Issue 8 - March 2022 Part No: MKT506\_USA  
Bedfont® Scientific Limited reserve the right to change or update this literature without prior notice.  
Registered in: England and Wales. Registered No: 1289798

**USA Representative for Sales, Service & Support:**



Official Bedfont Distributor  
212 Cottage Grove Ave Suite C  
Santa Barbara, CA 93101  
Tel: 800.707.5751  
Fax: 800.721.2377  
Email: [service@covita.net](mailto:service@covita.net)  
Website: [www.covita.net](http://www.covita.net)



**Stephen Rowe**  
Cristimar E4-1  
Ave Juan Carlos I  
Los Cristianos, Arona, 38650  
Santa Cruz de Tenerife, Spain



MD 502905