
Smokerlyzer® CO Chart



COppm	%СОНЬ
30	5.43
29	5.27
28	5.11
27	4.95
26	4.79
25	4.63
24	4.47
23	4.31
22	4.15
21	3.99
20	3.83
19	3.67
18	3.51
17	3.35
16	3.19
15	3.03
14	2.87
13	2.71
12	2.55
11	2.39
10	2.23
09	2.07
08	1.91
07	1.75
06	1.59
05	1.43
04	1.27
03	1.11
02	0.95
01	0.79

Reference:

1. COppm-%COHb calculation taken from: Jarvis M et al (1986) "low cost Carbon Monoxide monitors in smoking assessment." Thorax 41 pp 886-887. Issue 4- Novemberl 2015, Part No: LAB469