

GastroCHECK[®] Gastrolyzer[®] Maintenance Log

Calibration of the monitors should be performed <u>monthly</u> regardless of use or when prompted by the monitor.

Reminder: Schedule your 2 year service and 5 year service.

Scheduled Date	Actual Date	Calibration	Cleaning	Name Signature



www.covita.net | tel 800.707.5751 | fax 800.721.2377 212 Cottage Grove Ave. Suite C Santa Barbara, CA 93101

© Copyright 2023 coVita LLC All Rights Reserved v051623| Reproduction permitted for organization use. Log is for informational use only. coVita makes no claims regarding the accuracy of the content of this form or it's use in conjunction with any diagnosis.