



The NEW Micro⁺™ basic

Helping people to stop smoking, one breath at a time.




CO monitoring
made easy

FREE
COdata⁺
™
software

Ideal for in-depth stop
smoking advice delivered by:

- Stop smoking services
- General practice
- Schools
- Pharmacies
- Substance Abuse
- Wellness
- Researchers



Smokerlyzer[®]
www.covita.net | 800-707-5751

Features and Benefits

- *Save time - instant boot up*
- *Color Touch Screen - simple & intuitive user interface*
- *Automatically saves last 150 readings*
- *Higher range of 0-150 ppm*
- *Easier to use larger color touch screen*
- *SteriTouch® technology to kill harmful bacteria*
- *Brand new, easy to use simple interface*
- *Even simpler calibration using on-device step by step instructions*
- *Automatic firmware upgrades*

Technical Specification

Concentration range	0-150ppm
Display	Full color touchscreen
Detection principle	Electrochemical sensor
Repeatability	<±5%
Accuracy	2ppm/±5%*
Power	3 x AA (LR6 or equivalent) - up to 1000 minutes 1 x CR2032 Lithium coin cell
T ₉₀ response time	<30 seconds
Operating temperature	15-40°C
Storage/transport temperature	0-50°C
Operating/storage/transport pressure	Atmospheric ±10%
Operating humidity	15-90% non-condensing
Storage/transport humidity	0-95%
Sensor operating life	5 years
Sensor sensitivity	1ppm
Sensor drift	<5% per annum
Dimensions	Approx. 37 x 77 x 140mm
Weight	Approx. 215g (including batteries)
Materials	Case: polycarbonate/ABS blend with SteriTouch® anti-microbial additive D-piece™: polypropylene SteriBreath™: polypropylene
H ₂ cross interference	≤6%

*Whichever is greater



British Technology

coVita | Bedfont Scientific Ltd

212 Cottage Grove Ave, Suite CSanta
Barbara, CA 93101
Tel: (800) 707-5751 Fax: (800) 721-2377
Email: info@covita.net
www.covita.net

© Bedfont Scientific Limited 2017

Issue 7 - April 2017, Part No: MKT377
Bedfont Scientific Limited reserve the right to change or update this literature without prior notice.
Registered in: England and Wales. Registered No: 1289798

